

DISPENSER TABLETOP FOAM M8280



Touch
Free



Infrared
Sensor



Battery
Operated



Perfect
Dosage



Easy To
Refill



Place
Anywhere



0.25s
Start Up

PRODUCT INFORMATION

- + 280 ml / 9.5 fl Oz
- + 200 doses
- + 3 x AA batteries
- + Infrared sensor
- + Perfect dosage
- + Easy to refill
- + Place on any flat surface
- + 0.25 second start up time

INSTRUCTIONS

To be used with a non-alcoholic antibacterial liquid.



Place one hand under the dispenser.



The foam comes out instantly.



Apply the foam thoroughly on your hands.



Keep rubbing until the foam has dried in.

SAFETY



- + Unclassified formula
- + Non-flammable, non-explosive
- + Non-hazardous material and liquid
- + No risk for abuse

